

Alkaline Resistant Syringe



Suitable for alkaline solution

Chemical
Resistance

Strong alkaline resistance enables maintenance of precision of I.D. and long life

High pressure
Resistance

100,000 times piston motion. Plunger tip is replaceable

■ Feature

- Strong alkaline resistance enables maintenance of precision of I.D. and long life
- **Ceramic** syringe and **PEEK** syringe are available
- Volume accuracy and repeatability : within $\pm 1\%$
- Adequate mechanical strength avoids easy breakage by outer shock
- We design special shapes for tip of syringe and end of plunger
- Offer the optimal syringe with careful consideration of purpose and condition of use by collaborating with you from the design phase
- Order from 1 piece trial is available

◆ Product specifications are subject to change without prior notice for improvement.



Contact Japan: 3800-16, Obuchi, Fuji, Shizuoka, 417-0801 JAPAN
Email: info@exmire.com URL: <http://www.ito-ex.co.jp/en/>
Contact China: 中国北京市朝阳区霄云路32号2号楼707室
Email: chai@exmire.com URL: <http://www.ito-ex.co.jp/zh/>

