

# Special Auto Sampler Syringe



Syringe designed for various customer's instrument

High Quality

Special designed syringe based on reliable Ito Microsyringe

Various Specification

Wide variation with our processing technologies

High Endurance

100,000 times piston motion. Plunger tip is replaceable


## ■ Feature

- Ito Microsyringes are highly esteemed and used worldwide in the field of chromatograph, biochemistry, medical inspection and industrial instrument which require precise injection.
- Inner surface of annealed hard glass is mirror-like polished and very smooth. It enables stable precise injection and does not vary across the age.  
Accuracy and reproducibility : within  $\pm 1\%$
- Best specification suited to each customer's purpose of use and condition is able to be offered.  
We will cooperate with customer from the design phase.

Specification example (Head shape)



Specification example (Plunger head)



Please contact us for other special shapes.

◆ Product specifications are subject to change without prior notice for improvement.



Contact Japan: 3800-16, Obuchi, Fuji, Shizuoka, 417-0801 JAPAN  
Email: [info@exmire.com](mailto:info@exmire.com) URL: <http://www.ito-ex.co.jp/en/>  
Contact China: 中国北京市朝阳区霄云路32号2号楼707室  
Email: [chai@exmire.com](mailto:chai@exmire.com) URL: <http://www.ito-ex.co.jp/zh/>

